

Don't forget to schedule time in the LATIS/AISOS calendar!



[Follow along with the protocol as you watch the video](#)

## GETTING SETUP FOR GEOMAGIC ALIGNING


1. Navigate to the folder "NMNH Summer 2023 david scans (raw-not geomagic aligned-not fused)"
  - a. Choose whichever cabinet folder you are working in
  - b. Choose whichever box folder within the cabinet folder you are working in
  - c. Copy the folder name of the specimen you will be aligning
2. Navigate to the folder "NMNH Summer 2023 david scans (geomagic aligned-not fused)"
  - a. Create a new folder, paste the copied name from above into the folder and add \_cabinet#\_drawer#\_box#\_GA
    - i. Ex: SWT027-6\_cabinet1.b\_drawer1\_box2\_GA
3. Open the google drive excel file for the data: [NMNH Data Collection Summer 2023](#)
  - a. Complete the fields in columns A-F

## GEOMAGIC


1. Open Geomagic Design X
2. Either open the folder with the raw scan files and drag in all of the individual scans
  - a. OR **import** the scans by selecting import (top left hand of the screen, red arrow pointing down) and navigating to the proper raw scan file folder and using the CTRL+A to select all scan files.
  - b. **Change the sample size unit to mm for every scan if geomagic asks!**
3. On the left hand navigation panel, the bottom section labeled "Model", there will be an option labeled "**Mesher**" once the individual scan data has uploaded.
  - a. Click "**Mesher**" to automatically select all the meshes at the same time



4. Either in the “HOME” or “ALIGNMENT” tabs, click the icon labeled “Align Between Scan Data”





- Select the open circle next to “Global And Fine” and ensure the selections match the photo below
- Click the ✓ next to “METHOD”
- The value next to “Max. Average Deviation” (should be a 0.00# mm value), record that value in golem G of the spreadsheet. “alignment error (mm)”
- Click the ✓ next to “Align Between Scan Data” to complete the alignment
- If you forget to record the Max Average Deviation value before completing the alignment, that's OK! Just follow steps 3-4 again.



5. Make sure all the meshes are highlighted still! (if not, repeat step 3.a), click **export**

- Click the Check mark next to export



- Navigate to the folder in “NMNH Summer 2023 david scans (geomagic aligned-not fused)” you created for the aligned files in [“Getting Setup”](#)

- Click the box next to “Export Each Scan Data As A Separate File”
- Change the file name to **SWT###-##\_GA**
- Save the files as **.obj** files

6. Repeat these steps for every set of scans.



- If there's a scan file with a single scan, put N/A in columns F and G and Save the scan in the new folder following 5.b.ii
- If the scan looks funky after aligning, contact project overlord (Samantha).

## GETTING SETUP FOR DAVID FUSING


- Navigate to the folder “NMNH Summer 2023 david scans (geomagic aligned-not fused)”
  - Copy the folder name of whichever specimen you are working on
- Navigate to the folder “NMNH Summer 2023 david scans (FUSED)”
  - Create a new folder, paste the copied name from above into the folder and add **\_FUSED**
    - Ex: SWT027-6\_cabinet1.b\_drawer1\_box2\_FUSED




6. Open the david software
7. In the Shape Fusion tab (should automatically open to this tab), on the right hand side of the software, click the **+** within the box labeled **"List of Scans"**
  - a. Navigate to the folder "NMNH Summer 2023 david scans (geomagic aligned-not fused)" and open the folder for the specimen you are working with, Select all scans (Ctrl + a) and click "Open"
8. On the left hand side locate the box labeled **"Fusion"**
  - a. Use the +/- buttons to change the resolution so that the vertex spacing is as close to 0.05mm as possible
    - i. If the options are 0.6 and/or 0.4, go with 0.4.
  - b. Sharpness should be **0**
  - c. Close Holes box should be **UNCHECKED**.
  - d. Click fuse!
9. Open the google drive excel file for the data: [NMNH Data Collection Summer 2023](#)
  - a. Locate the specimen you are working with and complete columns H-L



10. On the right hand side locate the box labeled **"Fusion Results"** and **click save** (bottom right corner)
  - a. Navigate to the folder "NMNH Summer 2023 david scans (FUSED)"
    - i. Locate and open the folder you created in step 5
    - ii. Change the name from FusionResult to specimen ID\_fused
      1. Ex: SWT###\_##\_fused



11. On the right hand side locate the box labeled **"Fusion Results"** and click the **red X** (bottom left corner) to return to the individual scans
12. On the left hand side locate the box labeled **"Alignment"**, use the dropdown arrow to select Global Fine Registration
  - a. Set iterations to 80
  - b. Click "Align Scans"



13. Click fuse as you did in step 8
14. On the right hand side locate the box labeled **"Fusion Results"** and click save (bottom right corner)
  - a. Navigate to the folder "NMNH Summer 2023 david scans (FUSED)"
    - i. Locate and open the folder you created in step 5
    - ii. Change the name from FusionResult to specimen ID\_GFA\_fused
      1. Ex: SWT###\_##\_GFA\_fused